

# Dentatus

## PRODUCT PRICE LIST 2014

DENTATUS - A WORLD LEADER IN THE DESIGN AND MANUFACTURER OF DENTAL PRODUCTS, SINCE 1930.

*Whatever your crown retention requirements...*

**SURTEX  
UPGRADE 2014  
WITH NEW  
PACKAGING**

### DENTATUS CLASSIC SURTEX® POSTS HAS THE SOLUTION:

- Premium dental posts with proven reliability
- Surface treatment for superior retention
- Innovative and practical assortment box for ease-of-use
- Swedish high quality posts with 100% quality control



TITANIUM



GOLD PLATED



STAINLESS STEEL



MADE IN SWEDEN



POST SYSTEMS

---

OVERVIEW - POST SYSTEMS & CORRESPONDING REAMERS	4-5
DENTATUS CLASSIC SURTEX® POSTS	6-13
ACCESSORIES FOR DENTATUS METAL POSTS	13
LUSCENT ANCHORS®	14-16
TWIN LUSCENT ANCHORS®	18-19
LUMINEX® LIGHT TRANSMITTING POST	20-21

REAMERS

---

OVERVIEW - REAMERS	22
HELIX REAMERS (CLASSIC/CONICAL)	23-24
STANDARD REAMERS (CLASSIC/CONICAL)	25
SPECIAL REAMERS	26

FINISHING & POLISHING

---

PROFIN® DIRECTIONAL SYSTEM	27-29
PROFIN® BY HAND	29
LAMINEER TIPS FOR PROFIN/EVA	30-33
PER-IO-TOR ROOT PLANING INSTRUMENTS	34-35

MISCELLANEOUS

---

MATRIX BANDS	36-38
AMALGAM VIBRATOR PLUGGERS	39

ARTICULATORS & ACCESSORIES

---

MAAXLOC & ARTICULATORS	40-41
ARTICULATOR ACCESSORIES & MOUNTING PLATES	42
FACE BOWS & ACCESSORIES	43

NARROW BODY IMPLANT SYSTEMS

---

ATLAS	44
ANEW	45

GENERAL TERMS

---



## **LUSCENT ANCHORS®**

Luscent Anchors® non-metallic posts transmit light for curing composites within the deep confines of canals providing anchorage and flawless aesthetics. Resin reinforced fiberglass Luscent Anchors® fulfill dentistry's preference for very strong, non-metallic posts that are consistent with conventional adhesive techniques and materials to support core foundations and aesthetic restorations. The radioluscent anchors transmit tooth color without being visible under the most adverse light conditions. They also provide a sound core-crown gingival interface radiating natural color through soft and hard tissues surrounding teeth.

Corresponding reamers to your LUSCENT Anchors®

HELIX CONICAL: RSG #3-6, RSH #3-6, RSI #3-5

STANDARD CONICAL: RUG #1-6, RUH #4-6.



## **TWIN LUSCENT ANCHORS®**

Twin Luscent Anchors® offer additional advantages to those of Luscent Anchors® – improved mechanical retention, both against pulling and rotational forces. It gives a possibility to choose between tapered or parallel post shape, depending on which end is inserted in apical direction.

Corresponding reamers to your TWIN LUSCENT Anchors®

HELIX CLASSIC: RSA #3-6, RSB #3-6, RSC #3-5

STANDARD CLASSIC: RUA #3-6, RUB #3-6.



## **DENTATUS CLASSIC SURTEX® POSTS**

The Dentatus Surtex® Classic post has a micro-texturized metal finish, creating two-fold increase of retention compared to machined surfaces. The treated metal surface increases retention to acrylic resins and significantly reduces micro leakage. Their matte finish reduces glare for better aesthetics.

Corresponding reamers to your SURTEX® Classic posts:

HELIX CLASSIC: RSA #1-6, RSB #1-6, RSC #1-5

STANDARD CLASSIC: RUA #1-6, RUB #1-6.









## **LUMINEX® POST SYSTEM**

Luminex® light-transmitting plastic post polymerizes composites inside weak, debilitated root structures leaving an ideal space for subsequent placement of prefabricated or custom cast posts and cores. Either one of the Luminex Smooth and Grooved Plastic Posts can be used as burnout posts for custom cast posts and cores.

Corresponding reamers to your LUMINEX® posts:

HELIX CLASSIC: RSA #1-6, RSB #1-6, RSC #1-5

STANDARD CLASSIC: RUA #1-6, RUB #1-6.

Post system (post sizes)		Corresponding Reamers						
Surtex	Twin Luscent		<b>Helix</b>				<b>Standard</b>	
Aztec			Classic Reamers				Classic Reamers	
Luminex			Short 28 mm	Long 33 mm	X Long 37 mm		Short 28 mm	Long 33 mm
#1	-		RSA-1/3	RSB-1/3	RSC-1/2		RUA-1/3	RUB-1/3
#2	-		RSA-2/3	RSB-2/3	RSC-2/2		RUA-2/3	RUB-2/3
#3	White (XS)		RSA-3/3	RSB-3/3	RSC-3/2		RUA-3/3	RUB-3/3
#4	Yellow (S)		RSA-4/3	RSB-4/3	RSC-4/2		RUA-4/3	RUB-4/3
#5	Red (M)		RSA-5/3	RSB-5/3	RSC-5/2		RUA-5/3	RUB-5/3
#6	Blue (L)		RSA-6/3	RSB-6/3	-		RUA-6/3	RUB-6/3
Luscent			Conical Reamers				Conical Reamers	
Purple (XXXS)			-	-	-		RUG-1/2	
Green (XXS)			-	-	-		RUG-2/2	
White (XS)			RSG-3/3	RSH-3/3	RSI-3/2		RUG-3/2	
Yellow (S)			RSG-4/3	RSH-4/3	RSI-4/2		RUG-4/2	RUH-4/2
Red (M)			RSG-5/3	RSH-5/3	RSI-5/2		RUG-5/2	RUH-5/2
Blue (L)			RSG-6/3	RSH-6/3	-		RUG-6/2	RUH-6/2

# Dentatus CLASSIC SURTEX® POSTS

PREMIUM DENTAL POSTS WITH PROVEN RELIABILITY

*Whatever your crown retention requirements...*



## DENTATUS CLASSIC SURTEX® POSTS HAS THE SOLUTION:

- Surface treatment for superior retention
- Innovative and practical assortment box for ease-of-use
- Swedish high quality posts with 100% quality control

The original Dentatus Classic Surtex® Posts have become a global standard in dental anchors, and offer dentists a reliable and simple solution for post and core build-up. Our high-quality, pre-fabricated posts are dependable and economic for most endodontic treatment needs, available in a range of materials including pure titanium, stainless steel and gold-plated.

## PRIMARY BENEFITS

What ever your crown retention requirements, Dentatus Classic Surtex Post system has the solution. The system is thoroughly tested for retention of core materials in endodontically treated teeth. The posts are cemented permanently or temporarily in the prepared root canal and will offer improved retention and a non-glare appearance due to its texturized surface. The system consists of reamers and specially fitted posts in 6 different diameters and 5 different lengths.

The anatomically shaped reamers and posts correspond to the natural morphology of the root canal to minimize risk for root fractures.

## SUPERIOR MATERIALS FOR UNPARALLELED QUALITY


The Titanium posts are made of pure titanium and are recommended for most retentive purposes. Titanium is biologically highly compatible with surrounding tissue.

The Stainless Steel posts are recommended for all retentive purposes where extra high shear strength is required.

The Gold Plated posts are suitable for most retentive purposes, especially where an economical alternative is required (e.g. in temporary crowns).




**SURTEX  
UPGRADE 2014  
WITH NEW  
PACKAGING**



## RELIABLE, SAFE, AND EASY TO USE

Successful dental treatment is dependent on a number of prerequisites: a treatment plan, the right equipment and a deep understanding of the task at hand. At Dentatus we provide reliable, safe, high-quality products, which simplify the doctor's work and reduce the time the patient spends in the chair.

The Dentatus Classic Surtex® Post system is the result of the experience and knowledge gained as a trusted partner to the dental industry since 1932. It comprises a wide assortment of posts and reamers to retain posts in endodontically treated root canals and is suitable for most retentive work in root canals. A clear system of size markings on the reamers and assortment box makes the system easy to understand and re-ordering components simple. Once the optimal post size has been selected, the canal is successively enlarged to the correct diameter using the reamer system. The post is then simply cemented in place ready for the desired crown type.



## SUPERIOR RETENTION THROUGH DESIGN

Dentatus Classic Surtex Posts have been designed to deliver superior retention properties with all standard restoratives, including resin-based and amalgam materials. The classic threaded post increases retention, but our innovative Surtex surface texture delivers 200% more mechanical retention compared to non-treated posts\*. Introduced in 1998, this textured surface not only minimizes the risk of posts loosening, but also reduces light reflection, avoiding unsightly post exposure through translucent materials.

Additional mechanical retention can be obtained through the expandable crown section of the post, where up to four corner sections can be extended by use of the cross-cut key provided with the assortment kit (all post types except stainless steel).


\*Nergiz et al. The Journal of Prosthetic Dentistry, Vol 78, No 5, Nov 1997

## OUTSTANDING ASSORTMENT IN AN INNOVATIVE BOX

You will never need to think twice about selecting the correct post or reamer again, as Dentatus Classic Surtex Posts are housed in an innovative line assortment box - the One-stop Dispenser. It stores all the different posts, reamers and depth-marking rings by size, accessible via the transparent rotating box lid, which allows one section to be open at a time.

The lid also incorporates measuring devices for confirmation of length and diameter of posts and reamers.

The easily refilled assortment box gives you everything you need and when you need it in one reliable, time-saving system.




**NEW FEATURE!  
ANTI-SLIP MATERIAL  
ON THE BOTTOM,  
TO KEEP STEADY  
ON TABLE TOP**

## THE ORIGINAL ANATOMIC POST

The Dentatus Classic Surtex® Posts design, comprising a parallel shaft with tapered apical end, combines the advantages of traditional conical and tapered posts. The tapered apical end mirrors the anatomy of the tooth, avoiding the risk of dentinal micro fractures, while the parallel shaft optimizes retention. Threading on the post increases its surface area into the surrounding cement, anchoring it securely in place. Launched in 1932, the research results and experience gathered form a knowledge base unequalled in the industry ever since.

The anatomic design comes in a wide range of sizes, from 7.8 to 17 mm in length and between 1.05 to 1.8 mm in diameter, to fit teeth of any size. Several posts can be combined for teeth with multiple roots.



*Square crosscut head for core material adaptation and retention.*

*Parallel sides for passive insertion with thread pitch for secure retention to all cements.*

*Anatomical geometry for safe extra length of posts and for tapered roots.*

## STARTER KITS

Titanium			Sug. Retail Price €
Assortment Kit			190,30
TST-60	60	Assorted Posts (20 sizes, 3 posts of each)	
	6	Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)	
	1	Hollow Key	
	1	Cross Key	
	20	Stoppers	
Complete Starter Kit			257,00
TST-B	120	Assorted Posts (20 sizes, 6 posts of each)	
	6	Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)	
	1	Hollow Key	
	1	Cross Key	
	20	Stoppers	
Starter Refill Kit - no Reamers			223,00
TST	120	Assorted Posts (20 sizes, 6 posts of each)	
	1	Hollow Key	
	1	Cross Key	
	20	Stoppers	

# STARTER KITS

<i>Gold plated</i>		<i>Sug. Retail Price €</i>
<b>Assortment Kit</b>		130,20
RST-60	60 Assorted Posts (20 sizes, 3 posts of each) 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm) 1 Hollow Key 1 Cross Key 20 Stoppers	
<b>Complete Starter Kit</b>		199,50
RST-B	120 Assorted Posts (20 sizes, 6 posts of each) 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm) 1 Hollow Key 1 Cross Key 20 Stoppers	
<b>Starter Refill Kit - no Reamers</b>		135,20
RST	120 Assorted Posts (20 sizes, 6 posts of each) 1 Hollow Key 1 Cross Key 20 Stoppers	
<i>Stainless Steel</i>		<i>Sug. Retail Price €</i>
<b>Assortment Kit</b>		182,90
SST-60	60 Assorted Posts (20 sizes, 3 posts of each) 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm) 1 Hollow Key 1 Cross Key 20 Stoppers	
<b>Complete Starter Kit</b>		247,20
SST-B	120 Assorted Posts (20 sizes, 6 posts of each) 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm) 1 Hollow Key 1 Cross Key 20 Stoppers	
<b>Starter Refill Kit - no Reamers</b>		214,20
SST	120 Assorted Posts (20 sizes, 6 posts of each) 1 Hollow Key 1 Cross Key 20 Stoppers	

## REFILLS

<b>Titanium</b>		<i>Sug. Retail Price €</i>
<b>Short, 7.8 mm, 15 pcs/pkg</b>		24,50
TST-S1	Diameter 1.05 mm	
TST-S2	Diameter 1.20 mm	
TST-S3	Diameter 1.35 mm	
TST-S4	Diameter 1.50 mm	
TST-S5	Diameter 1.65 mm	
TST-S6	Diameter 1.80 mm	
<b>Medium, 9.3 mm, 15 pcs/pkg</b>		
TST-M1	Diameter 1.05 mm	
TST-M2	Diameter 1.20 mm	
TST-M3	Diameter 1.35 mm	
TST-M4	Diameter 1.50 mm	
TST-M5	Diameter 1.65 mm	
TST-M6	Diameter 1.80 mm	
<b>Long, 11.8 mm, 15 pcs/pkg</b>		
TST-L1	Diameter 1.05 mm	
TST-L2	Diameter 1.20 mm	
TST-L3	Diameter 1.35 mm	
TST-L4	Diameter 1.50 mm	
TST-L5	Diameter 1.65 mm	
TST-L6	Diameter 1.80 mm	
<b>X long, 14.2 mm, 15 pcs/pkg</b>		
TST-EL3	Diameter 1.35 mm	
TST-EL4	Diameter 1.50 mm	
TST-EL5	Diameter 1.65 mm	
TST-EL6	Diameter 1.80 mm	
<b>Super long, 17 mm, 15 pcs/pkg</b>		
TST-SL3	Diameter 1.35 mm	
TST-SL4	Diameter 1.50 mm	
TST-SL5	Diameter 1.65 mm	
TST-SL6	Diameter 1.80 mm	

# REFILLS

Gold plated		Sug. Retail Price €
Short, 7.8 mm, 15 pcs/pkg		13,90
RST-S1	Diameter 1.05 mm	
RST-S2	Diameter 1.20 mm	
RST-S3	Diameter 1.35 mm	
RST-S4	Diameter 1.50 mm	
RST-S5	Diameter 1.65 mm	
RST-S6	Diameter 1.80 mm	
Medium, 9.3 mm, 15 pcs/pkg		
RST-M1	Diameter 1.05 mm	
RST-M2	Diameter 1.20 mm	
RST-M3	Diameter 1.35 mm	
RST-M4	Diameter 1.50 mm	
RST-M5	Diameter 1.65 mm	
RST-M6	Diameter 1.80 mm	
Long, 11.8 mm, 15 pcs/pkg		
RST-L1	Diameter 1.05 mm	
RST-L2	Diameter 1.20 mm	
RST-L3	Diameter 1.35 mm	
RST-L4	Diameter 1.50 mm	
RST-L5	Diameter 1.65 mm	
RST-L6	Diameter 1.80 mm	
X long, 14.2 mm, 15 pcs/pkg		
RST-EL3	Diameter 1.35 mm	
RST-EL4	Diameter 1.50 mm	
RST-EL5	Diameter 1.65 mm	
RST-EL6	Diameter 1.80 mm	
Super long, 17 mm, 15 pcs/pkg		
RST-SL3	Diameter 1.35 mm	
RST-SL4	Diameter 1.50 mm	
RST-SL5	Diameter 1.65 mm	
RST-SL6	Diameter 1.80 mm	

REFILLS

Stainless steel		Sug. Retail Price €
Short, 7.8 mm, 12 pcs/pkg		23,50
SST-S1	Diameter 1.05 mm	
SST-S2	Diameter 1.20 mm	
SST-S3	Diameter 1.35 mm	
SST-S4	Diameter 1.50 mm	
SST-S5	Diameter 1.65 mm	
SST-S6	Diameter 1.80 mm	
Medium, 9.3 mm, 12 pcs/pkg		
SST-M1	Diameter 1.05 mm	
SST-M2	Diameter 1.20 mm	
SST-M3	Diameter 1.35 mm	
SST-M4	Diameter 1.50 mm	
SST-M5	Diameter 1.65 mm	
SST-M6	Diameter 1.80 mm	
Long, 11.8 mm, 12 pcs/pkg		
SST-L1	Diameter 1.05 mm	
SST-L2	Diameter 1.20 mm	
SST-L3	Diameter 1.35 mm	
SST-L4	Diameter 1.50 mm	
SST-L5	Diameter 1.65 mm	
SST-L6	Diameter 1.80 mm	
X long, 14.2 mm, 12 pcs/pkg		
SST-EL3	Diameter 1.35 mm	
SST-EL4	Diameter 1.50 mm	
SST-EL5	Diameter 1.65 mm	
SST-EL6	Diameter 1.80 mm	
Super long, 17 mm, 12 pcs/pkg		
SST-SL3	Diameter 1.35 mm	
SST-SL4	Diameter 1.50 mm	
SST-SL5	Diameter 1.65 mm	
SST-SL6	Diameter 1.80 mm	

ACCESSORIES FOR DENTATUS METAL POSTS

		Sug. Retail Price €
NLH	1 Hollow Key	8,10
NLK	1 Cross Key	10,40
DPT	Sure Lock Tweezers-Special Tweezers for the Classic/Surtex Post System	27,00



PREMIER PRODUCT FOR STRONG AESTHETIC RESTORATIONS

Dentatus Luscent Anchors offer the One-Step procedure for simultaneously curing of adhesives and composites within the confines of canals, creating a mono-bloc strength anchor-core support.

The resin-reinforced fiberglass Luscent Anchors are the optimal choice for anchorage and a cohesive foundation for flawless aesthetics. The anchors can be economically placed in the smallest and largest canals with off-the-shelf brand name adhesives and resins. Luscent Anchors transmit natural hues and tooth colors through the teeth and the surrounding gingival tissue.

AESTHETICS

The Dentatus Luscent Anchors transmit light into the confines of canals, polymerizing adhesives and composites attaching firmly the anchor with core foundation. The Luscent Anchors, with modern adhesive materials, overcome the disadvantages of metal post products especially in the highly aesthetic visible zones. It creates the ideal color background for crowns and the surrounding gingival tissue that is accomplished with conventional techniques and procedures.

ADHESION

The Luscent Anchors are composed of longitudinal fiberglass rods in a strong resin matrix with a modulus of elasticity well within the range of dentin. It provides the strength needed for supporting restorations. The anchor-core construction creates the ideal mono-bloc strength foundation for aesthetic restorations.

TECHNICAL INFORMATION

Luscent Anchors provide extremely high adhesion to composite materials. The glass fibres contain silica similar to the fillers of composite materials. The post and core is retentive and antirotational, conserving sound radicular and coronal structure.


Available in 6 sizes, XXXS-L

SIZES

XXXS LENGTH: 15 MM, Ø 1,0 MM		S LENGTH: 15 MM, Ø 1,45 MM*	
XXS LENGTH: 15 MM, Ø 1,15 MM		M LENGTH: 17 MM, Ø 1,6 MM*	
XS LENGTH: 15 MM, Ø 1,3 MM*		L LENGTH: 19 MM, Ø 1,75 MM*	

\*Helix Reamers available for these sizes

QUICK GUIDE FOR LUSCENT AND TWIN LUSCENT ANCHORS



## STARTER KITS

Sug. Retail  
Price €

### Starter kit 3, with Standard Conical Reamers

134,70

LUC-K3	15 Luscent Anchors, color coded depth markers (5 of each size: small, medium, large) 3 Assorted Standard Conical Reamers (one of each RUG-4, RUG-5, RUG-6) 1 Probos I Pathfinder #3 15 Core Forms-to-Fit
--------	---

## REFILLS

Sug. Retail  
Price €

### 5 pcs/pkg

33,30

LUC-L5	Luscent Anchors Large (blue depth markers) 6 assorted Forms-to-Fit
LUC-M5	Luscent Anchors Medium (red depth markers) 6 assorted Forms-to-Fit
LUC-S5	Luscent Anchors Small (yellow depth markers) 6 assorted Forms-to-Fit
LUC-XS5	Luscent Anchors Extra Small (white depth markers) 6 assorted Forms-to-Fit
LUC-XXS5	Luscent Anchors XX-Small (green depth markers) 6 assorted Forms-to-Fit
LUC-XXXS5	Luscent Anchors XXX-Small (purple depth markers) 6 assorted Forms-to-Fit

### 15 pcs/pkg

93,60

LUC-L15	Luscent Anchors Large (blue depth markers) 15 assorted Core Forms-to-Fit
LUC-M15	Luscent Anchors Medium (red depth markers) 15 assorted Core Forms-to-Fit
LUC-S15	Luscent Anchors Small (yellow depth markers) 15 assorted Core Forms-to-Fit
LUC-XS15	Luscent Anchors X-small (white depth markers) 15 assorted Core Forms-to-Fit
LUC-XXS15	Luscent Anchors XX-small (green depth markers) 15 assorted Core Forms-to-Fit
LUC-XXXS15	Luscent Anchors XXX-small (purple depth markers) 15 assorted Core Forms-to-Fit

### 30 pcs/pkg

175,60

LUC-L30	Luscent Anchors Large (blue depth marker) 30 assorted Core Forms-to-Fit
LUC-M30	Luscent Anchors Medium (red depth marker) 30 assorted Core Forms-to-Fit
LUC-S30	Luscent Anchors Small (yellow depth marker) 30 assorted Core Forms-to-Fit

## CORE FORMS-TO-FIT, REFILLS

Sug. Retail  
Price €

### 100 pcs/pkg

34,90

LX-F	Assorted
LX-FS	Small
LX-FM	Medium
LX-FL	Large

- Yes please!

I would like to receive the Dentatus Newsletter.

If you would like the next issue of our newsletter please inform us via e-mail [info@dentatus.se](mailto:info@dentatus.se) or telephone +46 8 546 509 00.

**Dentatus® Newsletter**  
MAY 2013

**Dear customer,**  
Your 1<sup>st</sup> edition of the Dentatus 2013 Newsletter is here!  
It has been an exciting start of 2013 for Dentatus with new product launches and a successful participation at this year's IDS exhibition in Cologne.

**Dentatus exhibited at the 35<sup>th</sup> International Dental Show in Cologne, March 11<sup>th</sup> - 16<sup>th</sup>**  
The biggest event in the dental business, IDS, took place in Cologne in the beginning of March and for Dentatus it was a great success. It was a great pleasure for our team to meet many of our current customers as well as new business acquaintances who took the time to visit our booth. We really appreciated the fruitful discussions about our products and business relationships.  
In case you missed us at IDS, please read more about the product launches below.

**The launch of Elypse Platforms**  
We were glad to introduce the expansion of the ANEW Narrow Body Implant System featuring the new Elypse Platforms. The new platform allows dentists to immediately retrofit patient's lower dentures with an option of later converting the restoration into a fixed prosthetic restoration - all while maintaining a comfortable and secure dental situation.  
The uniquely designed Elypse Platforms greatly enhances the product line by blending the best features of the proven Atlas Denture Comfort technology with the ANEW Implant System. The new feature was received with great enthusiasm by many of our customers, both new and established.

  
The Elypse Platforms are available as a kit including platforms in two different retention levels as well as necessary instruments to prepare the denture, and a dismantling driver. Corresponding refill packages are of course also available.  
The ANEW Implant System is an extraordinary complement to other implant systems. It enables practitioners to offer additional restorative solutions for patients with insufficient bone, limited interdental spaces or converging roots.  
Should you wish to launch Anew Narrow Body Implants or add the Elypse Platforms to your Anew range of products we will provide you with all the marketing materials you would need for a successful launch.  
*Anew Narrow Body Implants  
One Implant. More Options.*

**Dentatus continued the successful launch of Profin® PDX handpiece**  
The Profin system is the professional's choice for finishing, polishing and contouring of both natural and restored dental surfaces. The new Profin® PDX handpiece, which has tactile ergonomics, smooth glare free surface and a miniature head for best access into limiting areas, was launched to a bigger audience at the IDS exhibition. The new handpiece was very well perceived by many new customers and will soon be introduced in new markets.  
The PDX handpiece is available both as a single item and as a PDX/SK starter kit with a Profin PDX reciprocating handpiece and the five most commonly used Lamineey® tips for fine finishing and also 20 EVA Plastic tips for polishing.

  
Should you wish to launch the new Profin PDX handpiece in your market we will provide you with a marketing kit containing all you may need for success. Additionally, we now have new Profin promotion videos, available in different languages and a demo-kit for your sales representatives when demonstrating the Profin range to your customers.  
[www.dentatus.se/finishing-polishing/profin](http://www.dentatus.se/finishing-polishing/profin)

**New CE-Certificates are now available**  
This year Dentatus has changed its Notified Body and the new CE Certificates are now available. This means that products may be delivered with both 0470 and 0413 as the Notified Body code for some time. This is a logistic matter and it will not jeopardize the safety.



**Dentatus will move to new location in November 2013**  
We would like to bring to your attention that Dentatus will move the factory and head office to new and more modern facilities in November. The new site is situated not far from the current address outside Stockholm, Sweden. The address to the







### A WINNING COMBINATION OF LIGHT TRANSMISSION, ATTRACTIVE AESTHETICS AND TWICE THE RETENTION

Dentatus is raising the standard of excellence with the Twin Luscent Anchors. This innovative design is your visible assurance against accidental debonding of adhesive and resin-core materials. The slim mid-section creates a "physical choke." The vent groove eliminates air resin entrapment and prevents rotational dislocation. It all adds up to a winning combination of light transmission, attractive aesthetics and twice the retention.

#### FEATURES

#### BENEFITS

##### Light transmitting

Safe and easy handling. Now you can be sure that the root canal filling is completely polymerized!

##### Translucency

Aesthetics: Invisible inside the core. No dark shade.

##### Extra retention

Hour-glass shape provides improved retention of the post. The post cannot come out.

##### Solid structure

Anchor and Composite creates a very strong unit as materials function well together. With metal posts you have different types of materials; cement, metal, core-material.

##### Flexibility

Almost like dentine.  
Twin Luscent = Approx. 20 GPa (dentine = 18 GPa)  
Carbon = 120 GPa  
Titanium = 110 GPa

##### Fast restoration procedure

With Twin Luscent you work faster and more efficiently. The method is very simple.

##### No rotation of the post

Longitudinal groove provides anti-rotational interlocking. The posts will not rotate.

##### Tapered or parallel

Design of ends provides possibility to choose tapered or parallel apical post shape in root canal, depending on operators preference and individual root anatomy. The doctor has the opportunity to choose.

##### Colour coded

Fast size identification on tray with four different colours.

##### Quick Guide

Page 15

## STARTER KITS

Sug. Retail  
Price €

### Starter kit 1, with Standard Classic Reamers

146,40

LUT-K1	15 Twin Luscent Anchors, color coded depth markers (5 of each size: small, medium, large) 3 Assorted Standard Classic Reamers (one of each RUA-4, RUA-5 and RUB-6) 1 Probos I Pathfinder #4 1 Probos II reamer/router #4 15 Core Forms-to-Fit
--------	---

## REFILLS

Sug. Retail  
Price €

### 5 PCS/PKG

35,20

LUT-L5	Twin Luscent Anchors Large (blue depth markers) 6 assorted Forms-to-Fit
LUT-M5	Twin Luscent Anchors Medium (red depth markers) 6 assorted Forms-to-Fit
LUT-S5	Twin Luscent Anchors Small (yellow depth markers) 6 assorted Forms-to-Fit
LUT-XS5	Twin Luscent Anchors Extra Small (white depth markers) 6 assorted Forms-to-Fit

### 15 PCS/PKG

105,50

LUT-L15	Twin Luscent Anchors Large (blue depth markers) 15 assorted Core Forms-to-Fit
LUT-M15	Twin Luscent Anchors Medium (red depth markers) 15 assorted Core Forms-to-Fit
LUT-S15	Twin Luscent Anchors Small (yellow depth markers) 15 assorted Core Forms-to-Fit
LUT-XS15	Twin Luscent Anchors Extra Small (white depth markers) 15 assorted Core Forms-to-Fit

## SIZES

### Extra small

Colour: White  
Ø 1,26 mm  
Length 16,7 mm

### Small

Colour: Yellow  
Ø 1,40 mm  
Length 16,7 mm

### Medium

Colour: Red  
Ø 1,54 mm  
Length 16,7 mm

### Large

Colour: Blue  
Ø 1,68 mm  
Length 16,7 mm



# Dentatus LUMINEX® POST SYSTEM

THE PATENTED LIGHT TRANSMITTING POST TECHNOLOGY




## THE VISIBLE CURE FOR SAVING COMPROMISED ROOTS

Formerly, fragile, thin-walled and flared roots have presented major restorative problems due to fracturing and poor retention of posts. Often extraction has been the only alternative. Today there is an easy-to-use, effective solution to this problem.

Performed at chairside in one visit, the clear Luminex Light Transmitting Posts enable the polymerization of acid-etched, bonded, light-cured composites within the entire root canal to be used with a corresponding range of retentive posts.

## OVERVIEW OF THE CLINICAL PROCEDURE



### Canal preparation

Prepare the root canal to the desired size and depth with suitable Dentatus Standard Classic or Helix Classic reamers to fit a correspondingly sized Luminex smooth post.

### Centering and depth control

Insert the smooth Luminex post into the prepared canal to full depth and place a red stop for depth indication. Remove the luminex post. Acid-etch and bond according to the manufacturer's instructions for bonding system.

### Polymerization

Apply a low viscous light-curing composite resin into the canal. Place the luminex post to its full depth. Remove the excess material from the coronal region. Light-cure the post faciolingually according to the manufacturer's instructions for the composite.

### Restored full length canal

Remove the Luminex post by first twisting in place and then pulling in an upward motion. Before cementation, retentive, antirotational indents can be made in the canal with Probos II or similar instruments.

### Complete restoration

Select a suitable length post that corresponds in diameter with the Luminex post used. Verify its dimensions using the gauges. Apply the luting or bonding material to the post. Use the hollow key for passive post placement in the canal. After the cement/bonding material has fully set, the cross-key can be used to form metal post heads according to the desired finishing crown shape. For future canal access, the cross-key may be used to unscrew phosphate cemented prefabricated metal post. (The Luminex grooved posts are excellent to use as a pattern for custom cast post and cores that fit passively in the reconstructed canal.)

## LUMINEX ONE-STOP DISPENSER

- Reamers with ID grooves in sizes 1-6
- Elastic stoppers for length indication
- Measuring lengths in mm
- Confirming diameters of reamers and posts
- All posts and reamers are available in refills



## STARTER KITS

	KIT	Sug. Retail Price €
LST-T	Complete Assortment Kit 30 Luminex Smooth Posts - DIP (5 of each size: 1-6) 30 Luminex Grooved Posts - DCP (5 of each size: 1-6) 6 Assorted Standard Classic Reamers (1 RUA of each size), 1 Hollow Key, 1 Cross Key, 20 Stoppers	199,30

## REFILLS

		Sug. Retail Price €
	<b>20 Luminex Smooth Posts, (22 mm long)</b>	28,50
DIP-1	Diameter 1.05 mm	
DIP-2	Diameter 1.20 mm	
DIP-3	Diameter 1.35 mm	
DIP-4	Diameter 1.50 mm	
DIP-5	Diameter 1.65 mm	
DIP-6	Diameter 1.80 mm	
	<b>20 Luminex Grooved Burnout Posts, (22 mm long)</b>	47,00
DCP-1	Diameter 1.05 mm	
DCP-2	Diameter 1.20 mm	
DCP-3	Diameter 1.35 mm	
DCP-4	Diameter 1.50 mm	
DCP-5	Diameter 1.65 mm	
DCP-6	Diameter 1.80 mm	

# Dentatus REAMERS - OVERVIEW

---



## **HELIX REAMERS®**

---

Dentatus Helix Reamers have sharp cutting edges to create perfect seats for Dentatus root canal posts. The Helix reamers are available in both classic and conical shapes and in the same sizes, to be used instead of, or as alternatives to, the standard classic and conical reamers.

Made of toughened surgical quality stainless steel, the reamers are reliable and long-lasting in their performance, with high cutting capacity, smooth cutting action and improved transport of cutting debris.

Sizes are indicated by number of grooves on the shank. Number of grooves correspond to the same size number for post size on the Dentatus metal posts.

Use Helix Classic Reamers for Surtex, Twin Luscent, Aztec and Luminex posts, Helix Conical Reamers for Luscent.



## **STANDARD REAMERS®**

---

The Original Dentatus reamers. Made in the same high quality stainless steel as all our reamers, for long-lasting performance.

Sizes are indicated by number of grooves on the shank, which correspond to size number for post size on the Dentatus metal posts.

Use Standard Classic Reamers for Surtex, Twin Luscent, Aztec and Luminex posts, Standard Conical Reamers for Luscent.



## **SPECIAL REAMERS**


---

Dentatus Probos I Pathfinder is used for initial and safe removal of root filling material. Its square cross-cut working tip minimizes the risk for perforation of the dentinal walls. A built-in safety feature makes the Pathfinder to fracture next to the shank, if it gets stuck in the root canal.

Dentatus Probos II Reamer/Router has a special profile design which creates additional indents in the dentinal wall, when moved laterally. The indentations add extra surface to the dentinal wall, improving retention.


Sizes are indicated by number of grooves on the shank. Number of grooves correspond to the same size number for post size.

# HELIX CLASSIC REAMERS

ASSORTED PACKAGES		Sug. Retail Price €
		44,50
RSA-123	3 Helix Classic Reamers, Short (1 of each size 1, 2 ,3)	
RSA-345	3 Helix Classic Reamers, Short (1 of each size 3, 4 ,5)	
RSA-456	3 Helix Classic Reamers, Short (1 of each size 4, 5, 6)	
RSB-123	3 Helix Classic Reamers, Long (1 of each size 1, 2 ,3)	
RSB-345	3 Helix Classic Reamers, Long (1 of each size 3, 4 ,5)	
RSB-456	3 Helix Classic Reamers, Long (1 of each size 4, 5, 6)	
RSC-123	3 Helix Classic Reamers, X Long (1 of each size 1, 2 ,3)	
RSC-345	3 Helix Classic Reamers, X Long (1 of each size 3, 4 ,5)	

REFILLS			Sug. Retail Price €
Short, 28 mm, 3 pcs/pkg			44,50
RSA-1/3	XXX-Small (#1)	Diameter 1.05 mm	
RSA-2/3	XX-Small (#2)	Diameter 1.20 mm	
RSA-3/3	X-Small (#3)	Diameter 1.35 mm	
RSA-4/3	Small (#4)	Diameter 1.50 mm	
RSA-5/3	Medium (#5)	Diameter 1.65 mm	
RSA-6/3	Large (#6)	Diameter 1.80 mm	
Long, 33 mm, 3 pcs/pkg			
RSB-1/3	XXX-Small (#1)	Diameter 1.05 mm	
RSB-2/3	XX-Small (#2)	Diameter 1.20 mm	
RSB-3/3	X-Small (#3)	Diameter 1.35 mm	
RSB-4/3	Small (#4)	Diameter 1.50 mm	
RSB-5/3	Medium (#5)	Diameter 1.65 mm	
RSB-6/3	Large (#6)	Diameter 1.80 mm	
X Long, 37 mm, 2 pcs/pkg			29,90
RSC-1/2	XXX-Small (#1)	Diameter 1.05 mm	
RSC-2/2	XX-Small (#2)	Diameter 1.20 mm	
RSC-3/2	X-Small (#3)	Diameter 1.35 mm	
RSC-4/2	Small (#4)	Diameter 1.50 mm	
RSC-5/2	Medium (#5)	Diameter 1.65 mm	

# HELIX CONICAL REAMERS

ASSORTED PACKAGES		Sug. Retail Price €
RSG-345	3 Helix Conical Reamers, Short (1 of each size 3, 4 ,5)	 43,10
RSG-456	3 Helix Conical Reamers, Short (1 of each size 4, 5, 6)	
RSH-345	3 Helix Conical Reamers, Long (1 of each size 3, 4 ,5)	
RSH-456	3 Helix Conical Reamers, Long (1 of each size 4, 5, 6)	
RSI-345	3 Helix Conical Reamers, X Long (1 of each size 3, 4 ,5)	

REFILLS		
Short, 28 mm, 3 pcs/pkg		43,10
RSG-3/3	X-Small (#3)	
RSG-4/3	Small (#4)	
RSG-5/3	Medium (#5)	
RSG-6/3	Large (#6)	
Long, 33 mm, 3 pcs/pkg		
RSH-3/3	X-Small (#3)	
RSH-4/3	Small (#4)	
RSH-5/3	Medium (#5)	
RSH-6/3	Large (#6)	
X Long, 37 mm, 2 pcs/pkg		29,10
RSI-3/2	X-Small (#3)	
RSI-4/2	Small (#4)	
RSI-5/2	Medium (#5)	

# STANDARD CLASSIC REAMERS

ASSORTED PACKAGES		Sug. Retail Price €
RUA-1/6	6 Classic Reamers, Short (1 of each size)	51,30
RUB-1/6	6 Classic Reamers, Long (1 of each size)	



REFILLS			
Short, 28 mm, 3 pcs/pkg			27,00
RUA-1/3FP	XXX-Small (#1)	Diameter 1.05 mm	
RUA-2/3FP	XX-Small (#2)	Diameter 1.20 mm	
RUA-3/3FP	X-Small (#3)	Diameter 1.35 mm	
RUA-4/3FP	Small (#4)	Diameter 1.50 mm	
RUA-5/3FP	Medium (#5)	Diameter 1.65 mm	
RUA-6/3FP	Large (#6)	Diameter 1.80 mm	
Long, 33 mm, 3 pcs/pkg			
RUB-1/3FP	XXX-Small (#1)	Diameter 1.05 mm	
RUB-2/3FP	XX-Small (#2)	Diameter 1.20 mm	
RUB-3/3FP	X-Small (#3)	Diameter 1.35 mm	
RUB-4/3FP	Small (#4)	Diameter 1.50 mm	
RUB-5/3FP	Medium (#5)	Diameter 1.65 mm	
RUB-6/3FP	Large (#6)	Diameter 1.80 mm	
2 pcs/pkg, Depth marked 7, 9, 11 mm			36,90
RAM-4	Small (length 28 mm)	Diameter 1.50 mm	
RAM-5	Medium (length 28 mm)	Diameter 1.65 mm	
RBM-6	Large (length 33 mm)	Diameter 1.80 mm	

# STANDARD CONICAL REAMERS


		Sug. Retail Price €
Short, 28 mm, 2 pcs/pkg		27,00
RUG-1/2	XXX-Small (#1)	
RUG-2/2	XX-Small (#2)	
RUG-3/2	X-Small (#3)	
RUG-4/2	Small (#4)	
RUG-5/2	Medium (#5)	
RUG-6/2	Large (#6)	
Long, 33 mm, 2 pcs/pkg		
RUH-4/2	Small (#4)	
RUH-5/2	Medium (#5)	
RUH-6/2	Large (#6)	



# SPECIAL REAMERS

		Sug. Retail Price €
RUF-1/6	6 PROBOS® I, Pathfinder Assorted (1 of each size)	45,30
RUD	4 PROBOS® II, Reamer/Router Assorted (1 of each size)	44,00

## PROBOS I, Pathfinder

	33 mm, 3 pcs/pkg, refills	25,60
RUF-1/3	Diameter 0.60 mm (#1)	
RUF-2/3	Diameter 0.75 mm (#2)	
RUF-3/3	Diameter 0.90 mm (#3)	
RUF-4/3	Diameter 1.05 mm (#4)	
RUF-5/3	Diameter 1.20 mm (#5)	
RUF-6/3	Diameter 1.35 mm (#6)	

## PROBOS II, Reamer/Router

	28 mm, 3 pcs/pkg, refills	31,10
RUD-3/3	Diameter 1.35 mm (#3)	
RUD-4/3	Diameter 1.50 mm (#4)	
RUD-5/3	Diameter 1.65 mm (#5)	
RUD-6/3	Diameter 1.80 mm (#6)	

# MANDRELS

		Sug. Retail Price €
	12 pcs/pkg, refills ø 2,35 mm	41,90
MLR-RA	Dentatus Mandrels RA, length 24 mm	
MLR-HP	Dentatus Mandrels, HP, length 47 mm	



# **Dentatus** PROFIN® DIRECTIONAL SYSTEM

THE ULTIMATE INSTRUMENTATION FOR PROFESSIONAL FINISHING AND POLISHING

## **FOR PROFESSIONALS WORKING WITH:**

- Periodontics and preventive dentistry
- Prosthodontics
- Orthodontics
- Dental CAD/CAM

## **UNIQUE FEATURES WITH PROFIN:**

- Reciprocating movement, easy to control
- Free rotation 360° or locked-in position
- No heat build-up
- Provides superior surface finish on natural teeth and all restorations



**Profin®** PDX



Scan and read more about Profin

# Dentatus PROFIN® DIRECTIONAL SYSTEM

RECIPROCATING TOOLS FOR FINISHING AND POLISHING OF ALL TEETH AREAS

## ADVANTAGES WITH PROFIN PDX:

- Tactile ergonomics, good grip and handling properties
- Smooth glare free surface, low contrast colour
- Easy to clean
- Hard anodized surface, high scratch resistance
- Miniature head for best access into incapacious areas

The Profin Directional System, with precise and directionally controlled movement of a 1.2 mm reciprocating motion, is ideal for delicate and safe shaping, contouring, and polishing of the natural dentition as well as all types of restorations.

The handpiece and the dedicated and wide assortment of abrasive tips offer unlimited access to hard-to-reach intraoral areas, allowing reshaping of enamel and restorations made of porcelain, metal, composite and amalgam.

The Profin system is an excellent tool for refining interdental gingival and incisal embrasures to a natural form, contouring and individualizing anterior restorations, reshaping contacting tooth surfaces, and finenishing the margins of porcelain or metal restorations prior to cementation.

The tip assortment includes diamond coated Lamineer tips in various shapes and grits, Tungsten coated tips for gentle adjustment of margins between composite and enamel, and plastic tips for professional polishing and cleaning of any tooth surface.



## PROFIN STARTER KIT

Sug. Retail  
Price €

Contents/SK	Assorted Lamineer Diamond and Tungsten Plated Tips 1 Lamineer color coded Tip Stand Plastic and Wood Polishing Tip Samples
PDX/SK	Kit with contra-angle, for Standard E Electric or Air Motor 5 assorted Lamineer tips
PDWN/SK	Kit with contra-angle W&H WA-67/1,1A (with spray) 10 assorted Lamineer tips
PDWNL/SK	Kit with contra-angle W&H WA-67/1,1LT (with spray and light) 10 assorted Lamineer tips

530,60

953,60

1034,60



PDWNL

## PROFIN-ET STARTER KIT

Sug. Retail  
Price €

PT-1	Kit with 1 Profin-Et prophy-head, for Doriot fitting. Incl 1 LTA-S36, 1 LTA-S75, plastic and wooden tips samples
------	---

276,20



Profin-et

## CONTRA-ANGLE

Sug. Retail  
Price €

PDX	Contra-angle only, Standard E fitting for Electric or Air Motor
-----	---

374,50

## PROFIN® BY HAND

Designed specifically for minor excess material, this manual instrument complements the mechanical Profin contra-angle in those hard-to-reach areas, where extra delicate handling is requested.



	Intro Kit
LTA/HB	Profin Hand Instrument with 4 Lamineer Tips (LTA-50, LTA-S30, LTA-150, LTA-S100)

Sug. Retail  
Price €

126,20

Dentatus








LAMINEER® TIPS FOR PROFIN/EVA

RECIPROCATING TOOLS FOR FINISHING AND POLISHING OF ALL TEETH AREAS



The color-coded Lamineer® tips are designed to fit into normal, interdental spaces and accommodate the most extreme space requirements. These tips with various diamond particle grain sizes will remove gross excess materials as well as provide the finest finish and touch-ups. All tips are safe sided and can be sterilized using conventional methods.





LAMINEER TIPS

Cat#		
Shape		Recommended use
		Removal of gross overhangs and shaping of all restorative material - use tips with large particle size. For shaping and polishing all prosthetic materials and enamel. Final polishing - use tips with small particle size. Tungsten Carbide does not scuff enamel. For final finishing of composite filling margins.
		Same use as above. Smaller and thinner tip provide better access to narrow areas and closer to interdental contact points. Tungsten Carbide does not scuff enamel. For final finishing of composite filling margins.
		Coated on left or right side. For stripping and removal of contact points.
		For occlusal and concave surfaces.
		For finishing shoulders, bevels etc.
		For shaping diverging surfaces and incisal edges
		For opening incisal and gingival, approximal areas.

LAMINEER® TIPS

LAMINEER® KITS		Sug. Retail Price €
LTA/U	20 Lamineer Tips, Standard & Knife Edge, Assorted, 1 Lamineer Tip Stand with color code chart	469,60
LTA	10 Lamineer Tips, Standard and Knife Edge, Assorted, 1 Lamineer Tip Stand with color code chart	261,90
LTA/S	10 Lamineer Tips, Knife Edge, Assorted, 1 Lamineer Tip Stand with color coded chart	261,90

ACCESSORIES		Sug. Retail Price €
BLS/L	1 Lamineer Tip Stand with capacity for 20 tips	38,50
LTA-10	1 Cleaning stone for diamonds	4,30

LAMINEER® REFILLS AND ACCESSORIES		Not actual sizes on tips
REFILL PACKAGES (2 Tips per pack) Safe sided in alternate 180° position.		Sug. Retail Price €
Diamond, Standard Thin Spatula		45,60
LTA-15/2	15 Micron, violet	<p>For gross reduction of overhangs, reshaping and fine finishing. All restorative and luting materials.</p>   <p>For finising of composite filling margins.</p>
LTA-30/2	30 “ , white	
LTA-50/2	50 “ , yellow	
LTA-75/2	75 “ , gold	
LTA-100/2	100 “ , green	
LTA-150/2	150 “ , red	
LTA 39/2	Wolfram , Black Tungsten	
Diamond, knife edge short		45,60
LTA-S15/2	15 Micron, violet	<p>Shorter, narrower with knife edge for access to fine wedge shaped embrasures.</p>   <p>For finising of composite filling margins.</p>
LTA-S30/2	30 “ , white	
LTA-S50/2	50 “ , yellow	
LTA-S75/2	75 “ , gold	
LTA-S100/2	100 “ , green	
LTA-S150/2	150 “ , red	
LTA S36/2	Wolfram , Black Tungsten	

## DIAMOND, PADDLE

Sug. Retail  
Price €

For smoothing larger flat or convex surfaces

45,60

LTA-D30/2 30 Micron , white  
LTA-D50/2 50 " , yellow  
LTA-D100/2 100 " , green  
LTA-D150/2 150 " , red



## DIAMOND, SAW BLADE KNIFE EDGE

For opening incisal and gingival embrasures.

45,60

LTA-T30/2 30 Micron, white  
LTA-T50/2 50 " , yellow



## DIAMOND, ORTHO STRIPPING BLADES

45,60

LTA-L50/2 50 Micron, Left yellow  
LTA-R50/2 50 " Right yellow  
LTA-LR50 50 " Left + Right, pair  
LTA-L75/2 75 " Micron, Left gold  
LTA-R75/2 75 " Right gold  
LTA-LR75 75 " Left + Right, pair



## DIAMOND, SHOULDER/BEVEL

For reshaping precise inclines and angles of shoulders and bevels.

45,60

LTA-C50/2 50 Micron, yellow



## DIAMOND, CONCAVE

For extra fine protrusive and lateral adjustments on concave surfaces.

45,60

LTA-E50/2 50 Micron, yellow



POLISHING TIPS (PLASTIC AND WOODEN)


Sug. Retail  
Price €

41,30

EVA-123 -100 100 Wing-Shaped (Blue - large)



EVA-123/S-100 100 Wing-Shaped, (Pink - small)



36,90

EVA-2000/100 100 Universal Plastic Polishing Tips,  
Wedge shape, for cleaning of all  
approximal areas (green)



EVA-5000/100 100 Plastic Polishing Tips, Thin spatula,  
blue (Replaces EVA-5-100)



16,60

BL-6580 "Professional Mechanical Tooth Cleaning  
(PMTTC), Finishing and Minimally Invasive  
Treatment of Caries and Periodontal Di-  
seases – Materials, Methods and Effects".  
Written by Prof. Per Axelsson



# Dentatus PER-IO-TOR

A NEW CONCEPT FOR SCALING, ROOT-PLANING AND DEBRIDEMENT

## PREVENT FURTHER REMOVAL OF ROOT CEMENTUM ONCE SURFACE IS CLEAN

The PER-IO-TOR instruments are specially designed to optimize cleaning and planing of the rough root cementum and prevent further removal of root cementum once the surface is clean and smooth.

The secret is in the unique design of the PER-IO-TOR instruments, which have planing grooves set in a smooth surface.

### RECIPROCATING




### ROTARY




## PER-IO-TOR




-for Profin reciprocating handpiece




PER-IO-TOR 1 and 2 are designed for root grooves and furcation areas.



PER-IO-TOR 3 and 4 should be used on plane and convex approximal surfaces. No 3 is extra thin for very narrow interproximal areas, such as lower anteriors.



PER-IO-TOR 5 has a concave working surface for vertical use on convex buccal and lingual surfaces.




PER-IO-TOR 6 has a convex working surface for vertical or diagonal use on concave approximal surfaces.

## PER-IO-TOR


### ROTARY




-rotary instruments for handpiece (RA)



PER-IO-TOR 7R is a tapered cone for vertical use on root surfaces and especially in root grooves



PER-IO-TOR 8R is a parabolic cone for use on concave root cavities.



PER-IO-TOR 9R is a hyperbolic taper for furcation areas.

PERIODONTICS

FACTS

- All PER-IO-TOR instruments are made of high quality stainless steel and may be cleaned in ultrasonic baths, in thermal washers / disinfectors and in autoclaves at 135° C.
- PER-IO-TOR 1-4 may be used in all Profin and EVA handpieces.
- PER-IO-TOR 5-6 should be used in a Profin handpiece, since these instruments should be used in a locked position.
- PER-IO-TOR 7R-9R have a standard rotating right-angle fitting (RA).
- Degree of removal: approximately 0,005 mm at 2-3 minutes use until the root cementum is smooth.
- All PER-IO-TOR instruments should be used at approximately 10.000 rpm.
- The PER-IO-TOR instruments have been developed by Prof. Per Axelsson, Sweden, and are carefully documented. For references, please contact Dentatus.
- For more information on the Profin and EVA handpieces, please contact Dentatus.



PER-IO-TOR ROOT PLANING INSTRUMENTS

RECIPROCATING		Sug. Retail Price €
For use with the Profin Directional reciprocating contra-angle.		
KITS		
TOR-K	Assortment Pack, TOR #1 to #6	183,40
REFILLS		
TOR-1	1 PER-IO-TOR Tip, #1, Round with grooves (s)	30,60
TOR-2	1 PER-IO-TOR Tip, #2, Round with grooves (m)	30,60
TOR-3	1 PER-IO-TOR Tip, #3, Spatula with eyelet	30,60
TOR-4	1 PER-IO-TOR Tip, #4, Spatula with ridges	30,60
TOR-5	1 PER-IO-TOR Tip, #5, Concave	25,30
TOR-6	1 PER-IO-TOR Tip, #6, Convex	25,30
ROTARY		Sug. Retail Price €
For use with rotary contra-angles.		
REFILLS		
TOR-7R	1 PER-IO-TOR Tip, #7R, Tapered cone	27,80
TOR-8R	1 PER-IO-TOR Tip, #8R, Parabolic cone	
TOR-9R	1 PER-IO-TOR Tip, #9R, Hyperbolic cone	









**DENTATUS MATRIX BANDS ARE MANUFACTURED FROM THE HIGHEST QUALITY STAINLESS STEEL. THE STANDARD SHAPE BANDS ARE AVAILABLE IN TWO THICKNESSES: 0.050 AND 0.040 MM.**







The 0.040 bands are made of high quality stainless steel, and the thicker bands, 0.050, are made from stainless steel with special dead soft quality. The choice of steel quality gives the bands high strength and good adjustability. The bands provide very good prerequisites for well contoured and defined interproximal contacts. To reduce total interproximal matrix band thickness on mesio-occlusal and disto-occlusal restorations without reducing actual band thickness, the MoDo bands are specially cut to go under the contact on the intact side of the tooth. Dentatus quality has proven to be highly regarded by users all over the world for a long period of time.

## MATRIX BANDS - SHAPES

### Original Matrix bands



### Special Class V, ad modum Dr. Nyström



MATRIX BANDS


Stainless Dead Soft Steel 0,050 mm (0,002")		
12 bands in package		Sug. Retail Price €
MDD-1/DS	# 1	3,30
MDD-2/DS	# 2	
MDD-3/DS	# 3	
MDD-4/DS	# 4	
MDD-5/DS	# 5	
MDD-6/DS	# 6	
100 bands in package		Sug. Retail Price €
MDD-1/100FP	# 1	16,40
MDD-2/100FP	# 2	
MDD-3/100FP	# 3	
MDD-4/100FP	# 4	
MDD-5/100FP	# 5	
MDD-6/100FP	# 6	
Thin Stainless Steel 0,04 mm (0,0016")		
12 bands in package		Sug. Retail Price €
MDD-1RF/DS	# 1	3,30
MDD-2RF/DS	# 2	
MDD-3RF/DS	# 3	
MDD-4RF/DS	# 4	
MDD-5RF/DS	# 5	
MDD-6RF/DS	# 6	
Special Class V, ad modum Dr. Nyström 0,040 mm (0,0016") Stainless Steel		
12 bands in package		Sug. Retail Price €
MDD-N1/DS	N1	3,30
MDD-N3/DS	N3	
MDD-N4A/DS	N4A	
MDD-N4B/DS	N4B	
MDD-N5A/DS	N5A	
MDD-N5B/DS	N5B	

RETAINERS

Sug. Retail  
Price €

MEC	Universal, Straight, 1 pc	35,70
MEC-N1	Universal, Angled, Right, 1 pc	
MEC-N2	Universal, Angled, Left, 1 pc	
MEC-N	Universal, Angled, 2 pc, 1 each - Right, Left	69,80

*MED small size retainers have been discontinued.  
Contact our sales department.*



# AMALGAM VIBRATOR PLUGGERS

A wide variety of plugger shapes for the Dentatus Amalgam Vibrator (AVA) are still available.

AVA Contra-angle available from Sirona.



## PLUGGERS, SEPARATE

Sug. Retail  
Price €

AVA-1	Cylindrical Plugger, Diameter 1.2mm
AVA-2	Cylindrical Plugger, Diameter 1.5mm
AVA-3	Cylindrical Plugger, Diameter 2.1mm
AVA-4	Cylindrical Plugger, Diameter 2.7mm
AVA-5	Cylindrical Plugger, Diameter 3.5 mm
AVA-11	Heart Shaped Plugger, Diameter 4.6mm
AVA-12	Heart Shaped Plugger, Diameter 5.7mm
AVA-13	Heart Shaped Plugger, Diameter 6.6mm
AVA-14	Heart Shaped Plugger, Diameter 8.1mm
AVA-15	Heart Shaped Plugger, Diameter 10.1mm
AVA-21	Spatula Plugger, Diameter 3.0 mm
AVA-31	Spherical Plugger, Diameter 5.0 mm

20,70



## MAAXLOC

A revolutionary aid to identify and compensate for individual asymmetries between the right and left mandibular height before mounting of models in an articulator. Minimizes risks for unbalanced occlusal relations to being permanented in restorations, helps the dentist to reduce pain and reduces costs by rationalisation of treatments.



ARL

## FOR EFFICIENT AND ACCURATE STUDY OF OCCLUSAL RELATIONS

Dentatus Articulator type ARL is basically the same as type ARH but with another unique feature which is not found in any other articulator. The ARL Articulator has an adjustable Upper Jaw Member, and with the GAUGE BLOCK the articulator can be adjusted to accentuate mutual conformity of the upper and lower jaw members. This adjustment is necessary if the upper jaw member for some reason has changed its position and also to assure correct basic measurements when cooperating dentists and laboratories work with their own ARL articulators.



## ARTICULATORS & FACEBOWS

The Dentatus series of articulators and facebows have been designed to fulfill the requirements necessary for the efficient and accurate study of occlusal relations, diagnosis, treatment planning and for producing quality occlusal restorations in natural and artificial dentitions. Our articulators are crafted to the highest manufacturing standards for long lasting precision and function.



ARH

## A TRUE COMPANION


Dentatus Articulator ARH proves its value by being the most sold of all the Dentatus Articulator models ever. Being a precision instrument but yet with an unsurpassed sturdy and robust construction, the Dentatus ARH has become the standard in many countries globally. Materials used and surface finish makes our articulators stand up against everyday use for many years.

MAAXLOC		Sug. Retail Price €
MAX	Mandibular Axis Locator	1112,40

## ARTICULATORS

		Sug. Retail Price €
ARH	<b>Standard Type Semi-Adjustable Articulator, including</b> <ul style="list-style-type: none"> <li>1 Solid Incisal Table (ARH-G)</li> <li>1 Incisal Pin (ARH-22)</li> <li>1 Orbital Indicator (ARH-K)</li> <li>1 Long Support Pin (ARH-39)</li> <li>2 Zinc Mounting Plates (ARH-Z)</li> <li>2 Composite Reinforced Mounting Plates (ARH-ZC)</li> </ul>	797,70
ARL	<b>Standard Type Semi-Adjustable Articulator, including</b> <ul style="list-style-type: none"> <li>1 Solid Incisal Table (ARH-G)</li> <li>1 Incisal Pin (ARH-22)</li> <li>1 Orbital Indicator (ARH-K)</li> <li>1 Gauge Block (ARL-M)</li> <li>1 Long Support Pin (ARH-39)</li> <li>2 Zinc Mounting Plates (ARH-Z)</li> <li>2 Composite Reinforced Mounting Plates (ARH-ZC)</li> </ul>	887,60

PLANAS		Sug. Retail Price €
ARL/P	Articulator ad modum Planas	938,40
ARL/P-P	Plexiplate kit for Articulator ad modum Planas	71,70
AEP	Facebow ad modum Planas	321,80
AEP-3T	Bitefork kit (3 pcs) for Facebow ad modum Planas	137,30



## ARTICULATOR ACCESSORIES

		Sug. Retail Price €
ARH-167	1 Incisal Pin for ARH-P	41,10
ARH-175	1 Rubber Foot	1,60
ARH-22	1 Incisal Pin for ARH-G & ARH-EF	25,30
ARH-34	1 Lock Nut for Orbital Axis Plane Indicator	12,60
ARH-39	1 Long Support Rod	14,40
ARH-EF	1 Flexible Incisal Table	276,40
ARH-G	1 Solid Incisal Table	75,80
ARH-K	1 Orbital Axis Plane Indicator	25,30
ARH-P	1 Segmented/Adjustable Incisal Table	243,20
ARL-M	1 Gauge Block (For ARL)	99,60
ARS-18	1 Allen Wrench (For ARA, ARH, ARL, ARD)	8,10
BDA	1 Telescopic Bitefork Support	167,50
GDA	1 Mounting Cradle	60,00
MGA	1 Remounting Jig	139,00
SLH	1 Wooden Carrying Case for Articulators (All Models except ARA)	229,10

## MOUNTING PLATES

		Sug. Retail Price €
ARH-Z-2	Zinc Mounting Plate, 2 pcs	11,80
ARH-Z-28	Zinc Mounting Plate, 28 pcs	150,50
ARH-ZC-2	Reinforced Composite Mounting Plate, 2 pcs	9,60
ARH-ZC-28	Reinforced Composite Mounting Plate, 28 pcs	120,40



## FACE BOWS

		Sug. Retail Price €
AEB	<b>Facia Face Bow</b> <i>For use with all articulator Models, including</i> 1 Bitefork (AEB-G) 1 Orbital Pointer (AEB-18) 1 Clamp For Bitefork (AEB-Y) 1 Clamp For Orbital Pointer (AEB-R) 2 Standard Calibrated Condylar Rods (AEB-B)	321,80
AEE	<b>Earpiece Face Bow</b> <i>For use with all Articulator Models, including</i> 1 Bitefork (AEB-G) 1 Orbital Pointer (AEB-18) 1 Clamp for Bitefork (AEB-Y) 1 Clamp For Orbital Pointer (AEB-R) 2 Special Calibrated Condylar Rods with Blue Rounded Earpieces (AEE-E2)	321,80



AEB



AEE

## FACE BOW ACCESSORIES

		Sug. Retail Price €
AEB-18	1 Orbital Pointer	19,00
AEB-19	1 Lateral Support Rod	34,80
AEB-G	1 Bitefork (Standard)	44,30
AEB-J	1 Bitefork (Filled)	53,80
AEB-NO	2 Calibrated Condylar Rods with Extendible Pins	118,60
AEB-R	1 Locking Clamp for Orbital Pointer/or Lateral Support Rod	94,60
AEE-E2	2 Calibrated Condylar Rods with Blue Rounded Earpieces	53,80
AEB-B2	2 Calibrated Condylar Rod (standard)	53,80

## SERVICE!

All Dentatus articulators and facebows can be sent to our servicedepartment for repairs and refurbishment. Please contact us for shipping, handling and price.

*- The Simple Solution for Denture Comfort*


### AN IMPLANT SYSTEM THAT FITS!

- FITS THIN RIDGES
- FITS YOUR PATIENT'S BUDGET
- FITS YOUR PRACTICE


The narrow diameter ATLAS implant is designed to retain lower dentures for patients with moderate to severely resorbed bone. ATLAS Implants utilized with Tuf-Link<sup>®</sup> silicone reline material provide cushioned support and stabilization with gentle, firm retention while distributing the chewing force on the ridge and implants, all without housings, O-rings or adhesives.

The ATLAS minimally invasive procedure without a surgical flap dramatically reduces chair time and the simple technical steps make it vastly more affordable to many patients. Because the implants can be immediately loaded and cushioned, patients can walk into your office and leave within an hour or less with a comfortable, stable denture.

ATLAS one piece implants are available in 1.8, 2.2, 2.4 & 2.8 mm diameter and 7, 10, & 14 mm thread lengths.



### DENTURE COMFORT



### SEMINARS & WORKSHOPS

Dentatus arranges hands-on courses for the Atlas and Anew Implant Systems in different countries.

For more information about these courses or to sign-up, please contact Dentatus via

e-mail [info@dentatus.se](mailto:info@dentatus.se) or telephone +46-8-546 509 00.




## MINIMALLY INVASIVE NARROW BODY IMPLANTS FOR SCREW RETAINED PROSTHETICS

Anew narrow body implants advance the field of implantology with unique solutions. They are ideal for patients who have limited interdental spaces, insufficient bone level requiring otherwise lengthy bone augmentation procedures, or financial constraint.

Anew's prosthetic components provide patients with a fixed restoration at the time of placement. The non-hygroscopic screwcap abutment allows for implants to be easily monitored, altered and adapted to a permanent restoration.

The Anew's user-friendly components provide unprecedented opportunities for single, multiple, and large restorations, as they can be placed in narrow ridges, limiting, approximating roots and teeth, and in minimal occlusal spaces.

Anew is available in 1.8, 2.2, 2.4 diameters and 7, 10, & 14 mm thread lengths.



Anew resin screwcap assembled restorations, removable without force, prevent distortion and overcome metal screw fatigue and implant failures.




### MULTI-UNIT RESTORATIONS



*Photos courtesy Dr. Ziv Mazor*



### SINGLE TOOTH PROSTHETICS



*Photos courtesy Dr. Philip VanMeter*



### ELYPSE® FOR REMOVABLE DENTURE



*Photos courtesy Dr. Mark Iacobelli*

For more detailed information about our implant systems please visit our web site [www.dentatus.com](http://www.dentatus.com) or contact us via e-mail [info@dentatus.se](mailto:info@dentatus.se) or telephone +46-8-546 509 00.

# **GENERAL TERMS & CONDITIONS**

## **DENTATUS BUSINESS POLICY**

It is our policy to make every effort to genuinely satisfy our customers by providing consistently high quality products and services. Our products are designed to perform as claimed when used in the manner recommended in our technical instructions and manuals.

## **PRICES**

Prices are subject to change without notice. Possession of a price list does not constitute an offer to sell. All prices are in European currency, Euro €.

## **ORDERING**

When ordering products, please include complete catalog number.

## **MINIMUM ORDERS**

Orders totaling less than Euro 200.- will be subject to a special processing charge of Euro 30.-.

## **TERMS**

Against documents unless otherwise agreed upon. Past due balances are subject to a 2% per month interest. Delivered products remain under the property of Dentatus AB until complete payment has been made by the customer.

## **PRODUCT CHANGES**

Dentatus reserves the right to make product changes without prior notice. The customer is prohibited to remove or change Dentatus' trademarks and trade names, signs or any marks on any Dentatus product or packaging or to make any alternations in the construction or design of any product.

## **SHIPMENTS**

All shipments are free factory, and unless otherwise instructed, will be sent the most economical way.

## **CLAIMS**

Short claims must be made within 10 days from receipt of goods. Damaged goods must be reported to the carrier within 24 hours. Damage claim must be filled by consignee.

## **RETURNS**

Before goods is returned for credit or exchange, prior written approval including a RMA (Return to Manufacturer Authorization) number must be obtained.

Only standard merchandise of current manufacture and less than one year from date of original shipment will be accepted for credit. Merchandise returned in the original, unopened container and accompanied by the invoice number and date of purchase will be credited at the original purchase price. However, returned "opened" merchandise or merchandise requiring repackaging will be subject to a 15% restocking charge. Return of merchandise must be prepaid. Merchandise returned without prior written approval will be returned to the customer at the customer's expense.

## **WARRANTY**

Dentatus warrants its products to the original purchaser against defective workmanship and materials under normal procedure of use within the dental profession for a period of 90 days from date of sale. Warranty repair claims must have proof of purchase to be valid.

Dentatus AB, January 2014



